

CURRENT STATUS OF IMAGING MICROBIAL BIOFILMS IN THREE-DIMENSIONAL OPAQUE POROUS MEDIA USING X-RAY MICROTOMOGRAPHY

Dorthe Wildenschild¹, Gabriel Iltis¹, Ryan Armstrong¹, Yohan Davit², James Connolly³, Robin Gerlach³, and Brian Wood¹




¹ School of Chemical, Biological and Environmental Engineering
Oregon State University



² Institut de Mecanique des Fluides de Toulouse, Université de Toulouse, France



³ Center for Biofilm Engineering
Montana State University



Biofilms and Biofilm Architecture

- Soils and rocks – and engineering



- Bioclogging





- Bioremediation



- Microbial enhanced oil recovery

- We image:


- Biofilm surface morphology in a porous medium
- Not internal geometry and/or individual cells
- Spatial (and temporal) arrangement of biofilm in 3D porous media



*Images courtesy of the Center for Biofilm Engineering, Montana State University


Research Objectives

- To measure 3D architecture of biofilms in porous media (soft object embedded in a hard matrix)
 - Past approaches:
 - Destructive methods such as thin-sectioning
 - Two-dimensional micromodel systems (e.g. *Thiener et al., 2002*)
 - Numerous studies of 3D growth on flat substrate (no 3D porous medium) using CLSM – or on single grain only
 - Magnetic Resonance Microscopy/Nuclear Magnetic Resonance (e.g. *Seymour et al., 2004*) – limited resolution and long acquisition times
 - Nano-scale x-ray tomography observations of the biomass only (minus porous medium) (e.g. *Thieme et al., 2003*)
 - 3D studies using CLSM on index-matched media (Leis et al., 2006)
 - To provide 3D biofilm geometries for validation of existing theory and numerical models:
 - Ex.:
 - Individual-based Models (IbM), e.g. *Kreft et al., 1998, 2001*
 - IbM with lattice-Boltzmann model (*Graf von der Schulenburg et al., 2009*)



Why x-ray tomography?

- ❑ Simple answer: Three-dimensional information in opaque porous systems!
- ❑ Fast (~ 5 -10 mins per scan)
- ❑ High resolution ($\sim 5 \mu\text{m}$)
- ❑ No need to index-match
- ❑ “Non-destructive”
 - ❑ (to the porous medium)
- ❑ Potential to quantify:
 - ❑ structural arrangement
 - ❑ feedback with permeability/hydrodynamics
 - ❑ growth patterns/rates
 - ❑ mass transfer rates





Porter & Wildenschild (2009)

- flow conditions
- growth conditions
- microbial species


Phase contrast tomography

- Based on differences in refractive index (Snell's Law)
- Index-matching used with CLSM
- Can be used with x-rays as well
 - Momose et al. (1995, 1996) - medical imaging
- Work at Swiss Light Course (TOMCAT beam-line) and XRT Ltd. in Melbourne, Australia
- Relatively small variations in refractive index between water and biofilm



Horizontal slice of
biomass? in polystyrene
bead pack

Computed X-ray Microtomography (CMT)





Advanced Photon Source at Argonne National Lab
Advanced Light Source at Lawrence Berkeley Lab



Limitations of CMT for Imaging Biofilm



- X-ray cross-section of a biofilm is similar to that of water → need contrast agent
- Conventional contrast agents, e.g. potassium iodide, commercial medical contrast solutions (Fenestra, Isovue) diffuse readily into the biofilm
- X-ray exposure is expected to kill or severely inhibit biofilm growth
- Access to synchrotron x-ray sources (monochromatic radiation) is advantageous, but not unlimited
- Commercial micro CT systems do not support specific photoelectric edge-enhanced delineation



Diffuse interface

X-ray CMT Contrast Agents




- To date, two particle-based contrast agents have been used for imaging biofilms in porous media:
 - 1. Silver-coated hollow glass microspheres
 - *Ittis et al. (2011), Water Resources Research*
 - 2. Barium sulfate suspension
 - *Davit et al. (2010), Journal of Microscopy*



(Adrienne Phillips and James Connolly, MSU, 2010)

1. Silver Microsphere Contrast Agent


- Silver coated microspheres
 - ▣ $\sim 10\ \mu\text{m}$ diameter
 - ▣ Neutrally buoyant
 - ▣ High attenuation via Ag absorption edge at $\sim 26\ \text{keV}$
- Could use a variety of elemental contrast agents (w. photoelectric edges between $\sim 15\text{-}40\ \text{keV}$ such as Rb, Sr, Pd, Ag, I, Cs, Ba....)




Flow cell
column

Biofilm-
attached silver
microspheres


3D Biofilm Geometry



3D Biofilm Geometry




Rodriguez and Bishop (2007)




Shewanella oneidensis biofilm in glass bead pack (Iltis et al., 2011)


Segmentation and point-wrap algorithm (delineation of biofilm)



Biofilm with silver
microspheres attached in
glass bead pack
- light microscopy



Glass bead pack
(yellow) with biofilm-
associated silver
particles (blue)



(Avizo Fire[®])
PointWrap
representation of
biofilm

2D Biofilm Analysis – Ag particle approach

CMT

Light microscopy

Feature A



Feature B




	Surface Extraction		Point Wrap	
	Surface Area	Diff.	Surface Area	Diff.
	mm ²	%	mm ²	%
Feature A: LM	0.521	0	0.797	0
Feature A: CMT	0.495	-4.91	0.765	-4.13
Feature B: LM	0.639	0	0.999	0
Feature B: CMT	0.648	+1.01	1.017	+1.82

Additional details on our companion poster (tonight) – and in *Ittis et al. (2011)*

2. Barium Sulfate Contrast Agent

- Commercially available, medical-grade suspension
- $\sim 1\mu\text{m}$ sized particles
- Substitutes for aqueous phase while size-excluded from the biofilm phase
- Biofilm phase doped with KI
- *Davit et al. (2010)* used conventional CMT system
- Synchrotron:
 - Aqueous phase identified via barium absorption edge at $\sim 37.5\text{ keV}$
 - Biofilm phase via iodine edge at $\sim 33.5\text{ keV}$
 - Image subtraction – requires highly accurate registration






Davit et al. (2010)

Contrast Agent Comparison in 2D and 3D:

Gabe Iltis, Ryan Armstrong, Yohan Davit, and James Connelly

- Biofilm imaging involves “living” at the beam-line while the biofilm grows
- Compare microsphere and suspension approach in both 2D micromodels and 3D flow columns
- 2D: CMT vs. CLSM
- Data for Ag and BaSO₄ for same systems
- Results in enormous amounts of data - and image processing tedium..




Contrast Agent Comparison

Gray Scale Data


Barium sulfate

Silver-coated microspheres




Contrast Agent Comparison: Segmentation

Barium sulfate



Aqueous phase disappears


Silver-coated microspheres




Biofilm delineated by Ag particles

Contrast Agent 3-D Rendering


Barium sulfate



Coated microspheres



↑
Flow
Direction



Contrast Agent Advantages and Limitations

(Silver) Microspheres

□ Advantages:

- Microspheres (particles) can be added to the influent flow stream
- Surface attachment has occurred for all the biofilms tested so far

□ Disadvantages:

- Biofilm resolution is limited to the mean diameter of the particles
- Visualization/quantification is directly dependent on the quality of particle coverage on the biofilm

(Barium Sulfate) Suspension

□ Advantages:




- Suspension fills the hydraulically available pore space facilitating “easy” segmentation
- Good biofilm delineation
- Particle size is $\sim 1 \text{ } \mu\text{m}$

□ Disadvantages:

- High density and viscosity requires dilution and special care during addition to prevent dislodging of biofilm
- Potential compression of low-density biofilm?


Silver particles from CLSM
James Connolly, MSU 2010

Stained biofilm from CLSM
James Connolly, MSU 2010



Silver microspheres from CMT
overlain on rhodamine stained
biofilm from CLSM

Zoomed CMT/
CSLM image




CLSM image

Quantitative Studies

- Effect of flow rate on biofilm growth and porosity change
- Effect of bacterial species
- Measure pressure changes, oxygen concentrations etc.
- - see poster


Next!

- Combined quantification of biofilms and precipitate for SBR project on Precipitation at Solution-Solution Mixing Zones in Porous Media (Rick Colwell, PI)



Alpkvist et al. (2006)

- Numerical Modeling



Graf von der Schulenburg et al. (2009)

Conclusions

- The use of particle- and suspension-based x-ray tomography contrast agents facilitates quantitative imaging of biofilm in three-dimensional opaque porous media
- Architecture and spatial distribution can be obtained over many centimeters, for a variety of bacteria, and in a short time frame
- The technique needs refinement for application under a variety of conditions
- Implications for biofilm modeling:
 - Experimental data can be collected non-destructively for calibrating or validating models incorporating biofilm growth and related impacts on transport pathways/hydrodynamics
 - And: *“provide sufficient scientific understanding such that DOE sites would be able to incorporate coupled physical, chemical and biological processes into decision making for environmental remediation and long-term stewardship”*

Acknowledgments

- Mark Rivers

- ▣ Advanced Photon Source, Argonne National Laboratory, GeoSoilEnviroCARS (Sector 13)

- Alastair MacDowell and James Nasiatka

- ▣ Advanced Light Source, Lawrence Berkeley National Laboratory, Beam Line 8.3.2

- Mike Marsh

- ▣ Visualization Sciences Group, Houston, TX

- Gerald Debenest and Michel Quintard

- ▣ Institut de Mecanique des Fluides de Toulouse

- Danielle Jansik and Sassan Ostvar

- ▣ School of Chemical, Biological and Environmental Engineering, Oregon State University

- Funding

- ▣ The Department of Energy, Subsurface Biogeochemical Research Program), Grant No. DE-FG02-07ER64417 and Grant No. DE-FG02-09ER64734

